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Procedure-
Solid Samples Prepare the instrument as outlined in the operator's instruction manual. Determine the instrument blank. Log in a minimum of 3 Blank replicates. Press the Analyze button on the instrument's screen. After a short delay, the loading head slide-block will open. Place a 502-344 Nickel Basket into the loading head. Note: samples using automation should be placed in the appropriate autoloader position before starting the analysis sequence. Once the sequence has started, the automatic analysis will start and end automatically. 无锡进口分析仪供应商上海禹重实业有限公司致力于提供分析仪，竭诚为您服务。

Procedure–Powder/Chip Samples Determine the instrument blank. Press the Analyze button on the instrument screen, the lower electrode will close and the analysis sequence will start and end automatically. Repeat steps 2b through 2g a minimum of three times. Set the blank following the procedure outlined in the operator's instruction manual. Instrument calibration/drift correction. Log in a minimum of three Standard reps. Weigh approximately *** grams of a calibration/drift standard into a *** Nickel Capsule; enter the mass and standard identification into appropriate rep fields. If applicable, place the nickel capsule in the appropriate autoloader position. Note: LECO Reference Material *** do not require preparation. See preparation statement on the reference material certificate. Solid Standards may be used to calibrate when chip or powder standards are not available.

Sample Preparation Sampling and sample preparation is an important issue because traditional methods used to obtain samples for oxygen and nitrogen determination are different from those recommended for hydrogen, especially when sampling molten metal. The main difference in steel sampling procedures for oxygen/nitrogen and hydrogen is due to the mobility of hydrogen. Special precautions must be used when sampling for hydrogen. From molten steel and iron, a sample must be quickly quenched in cold water and chilled in a refrigerant such as liquefied nitrogen or a mixture of acetone and solid carbon dioxide in order to reduce losses of hydrogen from diffusion. Losses of oxygen and nitrogen from diffusion are not a problem. A sample that is taken for hydrogen and chilled in a refrigerant can also be used for oxygen and nitrogen determination. However, a sample that is typically taken for oxygen and/or nitrogen determination is not suitable for hydrogen determination due to hydrogen loss (diffusion). Surface contamination must be removed by filing or light grinding, using care not to overheat the sample.

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Procedure-

Solid Samples Press the Analyze button on the instrument screen again, the loading head slide-block will close and the lower electrode will open. Clean the upper and lower electrode manually, or, if applicable, remove the crucible and press the analyze button to clean with the automatic cleaner. Add approximately 0.05g of *** Graphite Powder to a *** Graphite Crucible. Firmly place the crucible on the lower electrode tip or appropriate autoloader position. Press the Analyze button on the instrument screen, the lower electrode will close and the analysis sequence will start and end automatically. 分析仪，就选上海禹重实业有限公司，让您满意，欢迎新老客户来电！无锡进口分析仪供应商

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